

COMPANION MICROPHONE

PRODUCT FEATURES

The Companion Microphone is an optional accessory for IntriCon's Lumen 155 hearing device — used to help listen to speech and other audio signals.

2.4 GHz Proprietary Audio Streaming Technology

Directional and Omni microphone capability

Volume Control

Rechargeable battery - with Micro USB charging port

LED light indicators

Easy usability for end users



SYSTEM SPECIFICATIONS

Battery	Rechargeable Lithium Polymer
Battery Life	Active: 9 hours; Inactive: 4 months
Operating Range	Typical range of 40 ft (12 meters)
Power Supply	Non-replaceable Lithium Polymer battery recharged with External wall adapter (Mfr PN: PSA03F-050Q1-R by Phihong Technology)
Wireless	2.4 GHz
Operating and Charging Temperature	32°F to 113°F (0°C to 45°C)
Storage Temperature	-4°F to 113°F (-20°C to 45°C)
Microphone	Switchable omni- and uni-directional

LED LIGHT BEHAVIOR

Function	LED Behavior
Power On — Press and hold power button. If released before 2 seconds will abort sequence.	
Turning Off — Press and hold power button. If released before 2 seconds will abort sequence.	
Charging	
Charged	
Low Battery	
Volume Up button pressed	
Volume Down button pressed	
Pairing button pressed	

For more information on IntriCon products, visit www.intricon.com
or e-mail hearinghealthsales@intricon.com

Contact Information:

INTRICON CORPORATION
1260 Red Fox Road
Saint Paul, MN 55112
Tel: 651.636.9770
Fax: 651.636.9503
Website: www.intricon.com
Email: hearinghealthsales@intricon.com

In Europe:
INTRICON GMBH
Kesselschmiedstrasse 10
D-85354 Freising, Germany
Telefon: 049-8161-4804-0
Telefax: 049-8161-4804-18
Email: service@intricon.de

